

EN TOOTH ACCIDENT - WHAT MUST BE DONE?

1. Stay calm and reassure your child.

If your child is **bleeding** from the mouth, have him/her **bite on one of the swabs** included in the kids emergency box and **cool the outside**, for example, with a cooling compress wrapped in a cloth (instant cold or cold hot/cold pack).

2. Locate the tooth / tooth fragment.

IMPORTANT: Do not touch the root surface and do not clean the tooth / tooth fragment.

3. Keep the tooth / tooth fragment moist immediately (the tooth / tooth fragment must not dry out).

OPTIMAL: Place the tooth / tooth fragment in a **professional tooth rescue box** (e.g. miradent SOS Zahnbox®). A professional tooth rescue box is proven to be the best storage for preserving the tissue (storage time up to 48 hours).

TEMPORARY SOLUTION: Alternatives such as cold, low-fat UHT milk (sterile milk) or physiological saline solution can only be used for a short time (storage time for a maximum of 30 minutes).

NEVER: Never keep a tooth / tooth fragment **moist with normal water** and **never** store it in the **oral cavity**.

DR. TILL'S ACTION TIP WITH THIS SET:



- Place the tooth / tooth fragment in the included **transport bag**.
- **Moisten with 2 – 3 drops** of the included **physiological saline solution**, **wrap up** carefully (squeezing the air out), and then **seal the pressure cap**.
- Keep for **30 minutes maximum**.

4. Visit a dental practice/clinic as quickly as possible.

Given the right conduct, **permanent teeth** can be very often saved. If a **milk tooth** breaks off or is knocked out, let the dental practice/clinic help you decide what to do next.

following the recommendations of:

Bundeszahnärztekammer, Bundesverband der Kinderzahnärzte

KNB_452015_Rev_2024-11-28